



# WHOLE BOVINE COLOSTRUM (FROZEN)

Bottle 1 Kg



## DESCRIPTION



The ECI powder of bovine colostrum is colostrum of first milking collected from European farms (free of Brucellosis, Tuberculosis, Leucosis). This product is 100% natural. Rich in maternal antibodies, it has a highly nutritious and healthy value.

ECI pays special care to the entire production process in order to ensure preservation of the biological activity of its products.

ANALYTICAL CONTENTS	
PROTEIN	14 %
FAT	6 %
CARBOHYDRATE	3 %
HUMIDITY	76 %
ENERGY VALUE (KCAL/100G)	120 Kcal
SPECIFIC CHARACTERISTIC	
IMMUNOGLOBULIN (Ig)	≥ 60 g/kg

Possibility of certified **IBR free** products  
(no antibodies against BHV1 virus)

*Life*  
WITH A +

Possibility  
IBR free

High  
quality

Minimum  
60g  
Immuno-  
globulin



## DIRECTIONS FOR USE



The following steps have to be strictly respected to preserve optimum efficiency of the product:

- **DEFROSTING:** THE PRODUCT MUST BE LEFT TO THAW IN A BUCKET OF WARM WATER (35°C TO 45°C)

- Do not defreeze the product in a microwave
- Do not defreeze the product in a bucket of water over 50°C

- **ADMINISTERING:**

The product has to be given to the calf in 3 meals with a feeding bottle or an oesophageal tube. The colostrum has to be around 35°C when administered.

- **DOSAGE:** 3 LITRES OF ECI COLOSTRUM

**1<sup>st</sup> MEAL :**

- 1 litre of colostrum as quickly as possible after the calf's birth

**2<sup>nd</sup> MEAL :**

- 1 litre of colostrum 6 hours after the calf's birth

**3<sup>rd</sup> MEAL :**

- 1 litre of colostrum 12 hours after the calf's birth



### PACKAGING

1 kg Bottle

### PRESERVING

30 months at -18°C

48 hours at +4°C (fridge)

### CERTIFICATION/QUALITY INSURANCE



### CONTACT

**ECI S.A. – European Colostrum Industry**

Rue de la Fontaine, 17A • 6900 Marloie (Belgique)

Tel/Fax. : +32 84 41 31 50 • Email : [info@eci-colostrum.com](mailto:info@eci-colostrum.com)

Unit of Establishment : 2.204.883.442 • VAT N° BE 0841 637 326

**[www.eci-colostrum.be](http://www.eci-colostrum.be)**